



OTTOMED ENDOSCOPY SYSTEMS & ANCILLARIES


Better Technology, Better Results.







Branch Offices across INDIA



Pan-India Service Assist

Fully equipped Service Centres across the country.

Smart Team

Trained & Experienced
Service Engineers

Service Support

Quick & On-time response.
Online & On-location.

Training Programs

DNT Programs to keep you updated.

We are there for you..



About Us

Ottomed is a reputed and established name in the Healthcare Industry in India involved in manufacturing and marketing of healthcare devices with several product patents in our name. We are the only Indian manufacturers of flexible Endoscopes and related Ancillary products.

Our Vision: To enable Lifecare through sustainable Technologies.

Our Mission: To answer customer needs and practice high degree of corporate governance.

Our Culture: Trust, Transparency & Collaboration

We believe that constant innovation is the only way forward. Our state-of-art Research & Development (R&D) department designs and develops next generation products and devices with sustainable technologies based upon customer feedback and needs to enable healthcare professionals provide high quality lifecare.

Our strong focus in R & D has made us the world leaders in FIBRELESS ILLUMINATION LED-AT- TIP technology which puts an end to light-guide glass fibres in video endoscopes. Coupled with a focus on quality and speed to market, our aim to exceed market expectations for performance and value. Our manufacturing processes adhere to the highest global and GMP standards.

Our strong sales & service presence, with offices & service points across the country, offer quality and time bound services to all our customers.

Our high-quality products are available internationally in many countries.

Ottomed Endoscopy is a division of Mitra Medical Services LLP, a part of the Mitra Group.

INDEX

Sr.No	Content	Page No
1	About Us	01
2	Global Presence	03
3	Patents	04
4	Go Fibreless	05
5	Ottomed mBLU Technology	06
6	Tactile Insertion Tube	07
7	Auxilliary Jet	08
8	HD Video Processor SmartEye-II Plus	09
9	Video mBLU Gastroscope OVG-mB2H	11
10	Video mBLU Colonoscope OVC-mB2H	13
11	Video Paediatric mBLU Colonoscope OVC-PmB2H	15
12	Video Duodenoscope OVD-2H	17
13	HD Video Processor SmartEye-II	19
14	Video Gastroscope OVG-2H	21
15	Video Ultra-Slim Gastroscope OVG-2HUS	23
16	Video Colonoscope OVC-2H	25
17	Video Colonoscope OVC-P2H	27
18	Video Duodenoscope OVD-2H	29
19	HD Video Processor SmartEye-S3	31
20	Video Gastroscope OVG-S2+	33
21	Video Gastroscope OVG-S2+US	35
22	Video Colonoscope OVC-S2+	37
23	Video Colonoscope OVC-PS2+L	39
24	Video Bronchoscope OVB-S2+/ OVB Xs2+/ OVB-PS2+/ OVB-XPS2+/ OVB-CS2+	
25	Video Pleurascope OTS-S2+	43
26	Video Laryngoscope OVL-S2+	45
27	Video Ureteroscope OUR-S2/ OUR-S2+	47
28	Video Cystoscope OCY-S2+	49
29	HD Video Processor Insight II	51
30	Video Bronchoscope OVB-IN+/ OVB-PIN+/ OVB-XPIN+/ OVB-CIN+	53
31	Video Pleurascope OTS-IN+	55
32	Video Laryngoscope OVL-IN+	57
33	Video Ureteroscope OUR-IN+	59
34	Video Cystoscope OCY-IN+	61
35	Video Cholangioscope OBS-IN+	63
36	Ottomed Video Veterinary Endoscope System	65
37	HD Video Processor SmartEye-S3	67
38	Video Veterinary Endoscope OVV-S2+ (I/S)	69
39	Video Ultra-Slim Veterinary Endoscope OVV-S2+US (I/S)	71
40	HD Video Processor SmartEye-II	73
41	Video Veterinary Endoscope OVV-2H (L/I/S)	75 75
42	Full HD Laparoscopy Camera System	73 77
42	EndoPlus HD Video Camera OVP-EP	79
44	LED Light source OLS	80
77	LED LIGHT SOUTO OLO	00

INDEX

Sr.No	Content	Page No.
45	Ancillaries	81
46	Endocarbon-Illi	83
47	Endocarbon-III	85
48	Endopump II	87
49	Endowasher Qubey-III (OEW-3)	89
50	Endowasher Qubey-II (OEW-2)	91
51	Endowasher Qubey-T	93
52	Endohanger 2	95
53	Endohanger	96
54	EndoTrolley	97
55	Medical Grade Monitor-OVM-24HD-A	99
56	EndoClinic	101
57	EndoClinic - Medical Image & Records Management System	103-108
58	Non-Cardiac Metallic Stents	109
59	Blood Bags & Blood Banking	110
60	Peritoneal Dialysis Systems	111
61	Our Product lines and Certificates	112



Video Endoscopy Processors



Video Endoscopes







Endoscopy Accessories



EndoClinic





Global Presence Patents

- CIS: Kazakhstan, Kyrgyzstan, Russia, Tajikistan, Uzbekistan. Qatar, Israel, Dubai, Syria, tine, Iraq, Oman. Africa: Ethiopia, Kenya, Senegal, Malawi, Mali, Zambia, Mauritius, Burundi, South Africa, Liberia, Morocco, Mozambique, Nigeria, Djibouti, Botswana, Seychelles.
- **Guyana**, Latin America: Argentina, Brazil, Colombia, Guatemala.
 - Asia: Afghanistan, Sri Lanka, Nepal, Malaysia, Singapore, Thailand, Philippines, Bangladesh.



03





Optical Imaging System













GO Fiberless!

Introduction to MBLU Experience



With Ottomed Advanced Tech Fiberless LED-AT-TIP




Visible sharp detection, lively and crystal clear image



A new 1:1 torque sensitive tactile insertion tube (T-I-T) design.

Experience Enhanced Image Endoscopy





Dedicated Glass Light Guide Illumination bundle will always get damaged

Dark and Lifeless Image

Poor Detection and Diagnosis





White Light Image The Next-Generation Image Technology



mBLU1 Image Support Study Of Muscosal Color Characterization

Going Fiberless helps cut repair cost and downtime

LED vs XENON Light Source		
	LED	XENON
Working Life (approx)	5000 Hours	500 Hours
Energy Consumption	Low	High
Replacement Cost	Low Due to Longer Life	High Due to Shorter Life Span requiring more replacements

Good Design Provides A Great User Experience.






Tactile Insertion Tube (TIT)

Auxilliary Water Jet

A new 1:1 torque sensitive tactile insertion tube (T-I-T) design



The Quick Torque Transmission (QTT) technology in these scopes provide the endoscopists a rapid control tool to carry out minimally invasive, precise and effective procedures and enabling them to make truly-informed diagnostic and treatment decisions in the best interest of their patients.



07

A dedicated water irrigation channel to precisely wash the colonic mucosa area in the GI tract and gain better visibility for easier diagnosis and detection.







SmartEye-II Plus

Model: OVP-SEII Plus

Next Generation Digital HD Video Endoscopy System

SmartEye-II Plus is a Next Generation Digital HD Video Endoscopy System offering best-in-class features to Gastroenterologists. The HD platform makes diagnostic and therapeutic endoscopic treatments smooth & hurdle free.

09

*Compatible with all mB2H and 2H Series Endoscopes

Key Features









HD Imaging

mBLU

Extreme Close up Endoscopy (eCUE)

Multiple O/P & I/P Video Signal







PIP/POP

Touch Screen

Compact & Light Weight

Power Supply	Voltage	100-240 VAC		
Size	Voltage fluctuation	Within+10%		
Size	Frequency	50/60 Hz		
	. ,	Within+3 Hz		
	Frequency fluctuation Consumption electric power			
	Fuse rating	75VA		
	Fuse size	4 A, 250 V ø5.0 x 20.0 mm		
	Dimensions & Weight	95.0 x 20.0 mm 320(W) x 500(L) x 152(H) mm 10.3kg		
	Differsions & Weight	320(W) X 300(L) X 132(H) HIIII 10.3Kg		
Classification (medical electrical equipment)	Type of protection against electric shock. Degree of protection against electric stock of applied part. Degree of protection against explosion.	Class 1		
Observation	Video signal output	VBS composite (PAL), Y/C, HDMI, HD-SDI, DVI-D Simultaneous output is possible.		
observation	Video signal input	VBS composite (PAL), HDMI		
•	White balance	By pressing the white balance switch on the front panel, automatic white balance can be performed.		
	Color tone adjustment	The following color tone adjustments may be using the color tone adjustment on the "System Setup" menu * "R" control:±10 steps * "B" control:±10 steps.		
	Iris area selection	The area over which light is measured can be changed according the body part that is being observed. The following to choices can be made using the iris mode selector on the front pa Average: Normal observation Peak: When focusing on and/or observing a small bright are		
	Image Enhancement mode setting	Structure enhancement in the endoscopic images can be enhanced electrically to increase the image sharpness.		
	Gamma Correction	Gamma correction adjustments can be performed by adjusting the gamma value.		
	Switching the enhancement modes	The enhancement level can be selected from 3 levels enhancement modes (Low, Med. and High) using theimage enhancement mode on the front button panel.		
	Freeze screen display	The endoscopic image can be displayed as a frozen image, using the freeze switch key from endoscope.		
	Brightness Selector	The brightness can be adjusted using the brightness adjustment on the "system setup/camera adjustment" menu.		
	Remote control	The following function can be controlled by the endoscope's remote switches (Freeze, Record Light, White Balance, mBLU, Enhancement, Contrast, Focus, Iris, Image Size, Remove Data EXT In)		
	mBLU	This mode is used to switch the normal mode to mBLU mode with white or blue light or vice versa.		
	External Storage	USB Memory Stick		
Documentation	Patient Data	The following data and modes can be displayed on the video monitor using the keyboard-ID number,patient Name, Sex, Age, Date of Birth, & Comments.		
	Display the record state	Date/time (built in clock)		
Imago Storago	Pre-procedural patient data	For a maximum of 40 patients- ID number, patient Name, Sex, Age, & Date of Birth The following setting on the front panel are stored even when the video processor is turned		
Image Storage & Retrieval	Memorization of selected setting	OFF- System Setup, Setting the image record, Patient Data, Colour Tone, & White Balance.		
Memory Backup	Lithium Battery	Life: 5 years		
Illumination	LED (White)	At Endoscope Tip		
Air Feeding	Pump	Diaphragm type pump		
Method	Airflow	"Low", "High", "Stop"		
	Airflow Pressurization of water container (detachable)	Air pressurization of water container (detachable)		





HD VIDEO GASTROSCOPE

mBLU HD Video Endoscope



Blue Light Imaging



Model: OVG-mB2H (Compatible Video Processor: OVP-SEII Plus)

Key Features



Fiberless LED-At-Tip



Auxiliary Water Jet Channel



HD Imaging



Single Step, Waterproof Connector









Optical System	Field Of View	140°
Optical System	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
	Optical Technology - mBLU	mBLU Mode 1 & mBLU Mode 2
Insertion Tube	Distal End Outer Diameter	9.8 mm
	Insertion Tube Outer Diameter	9.5 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W/B)
		3. Objective Lens
		4. Instrument Channel Outlet
		5. Auxiliary Water Channel
Instrument Channel	Working Length	1055 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		30 cm 3/s
Bending Section	Angulation Range	UP 210°
		DOWN 90°
		RIGHT 100°
		LEFT 100°
Total Length		1365 mm
Compatible Video Processor		OVP-SEII Plus





HD VIDEO COLONOSCOPE



Model: OVC-mB2H (L/I/S) (Compatible Video Processor: OVP-SEII Plus)

Key Features



Fiberless LED-At-Tip



HD Imaging

13




Tactile Insertion Tube




Auxiliary Water Jet Channel



Single Step, Waterproof Connector









Specifications:


Optical System	Field Of View	140°
optical System	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
	Optical Technology - mBLU	mBLU Mode 1 & mBLU Mode 2
Insertion Tube	Distal End Outer Diameter	13 mm
	Insertion Tube Outer Diameter	12.8 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W/B)
		3. Objective Lens
		4. Instrument Channel Outlet
		5. Auxiliary Water Channel
	Working Length	L: 1680 mm
		I: 1330 mm
		S: 730 mm
Instrument Channel	Channel Inner Diameter	3.7 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		25 cm %
Bending Section	Angulation Range	UP 180°
		DOWN 180°
		RIGHT 160°
		LEFT 160°
Total Length		L: 1680 mm
		I: 1330 mm
		S: 730 mm
Compatible Video Processor		OVP-SEII Plus





PAEDIATRIC HD VIDEO COLONOSCOPE

mBLU HD Video Endoscope



Model: OVC-PmB2H (L/I) (Compatible Video Processor: OVP-SEII Plus)

Key Features










Fiberless LED-At-Tip

HD Imaging


Tactile Insertion Tube





Single Step, Waterproof Connector

000










Optical System	Field Of View	140°
	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
	Optical Technology - mBLU	mBLU Mode 1 & mBLU Mode 2
Insertion Tube	Distal End Outer Diameter	11.8 mm
	Insertion Tube Outer Diameter	11.8 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W/B)
		3. Objective Lens
		4. Instrument Channel Outlet
		5. Auxiliary Water Channel
	Working Length	L: 1680 mm
		I: 1330 mm
Instrument Channel	Channel Inner Diameter	3.2 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		25 cm %s
Bending Section	Angulation Range	UP 180°
		DOWN 180°
		RIGHT 160°
		LEFT 160°
Total Length		L: 1980 mm
		I: 1630 mm
Compatible Video Processor		OVP-SEII Plus





HD VIDEO DUODENOSCOPE

Excellent image quality & easy operability in Therapeutic procedure



Model: OVD-2H (Compatible Video Processor: OVP-SEII Plus, OVP-SEII)

17

Key Features



Fiberless LED-At-Tip



Single Step, Waterproof Connector









Specifications:

Optical System	Field Of View	100°
optical system	Direction Of View	Backward side viewing 5°
	Depth Of Field	5-60 mm
	Optical Technology - mBLU	mBLU Mode 1
Insertion Tube	Distal End Outer Diameter	13.5 mm
	Insertion Tube Outer Diameter	12.5 mm
	Distal End Enlarged	1. Air/Water Nozzle
	2 2 4	2. Objective Lens
	15	3. Illumination LEDs (W)
		4. Elevator
		5. Instrument Channel Outlet
Instrument Channel	Working Length	1260 mm
	Channel Inner Diameter	4.2 mm
	Minimum Visible Distance	10 mm from the distal end
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		25 cm 3/s
Bending Section	Angulation Range	UP 120°
		DOWN 90°
		RIGHT 110°
		LEFT 90°
Total Length		1545 mm
Compatible Video Processor		OVP-SEII Plus, OVP-SEII







SmartEye-II

Model: OVP-SEII

High Definition Video Endoscopy System

SmartEye-II is a High Definition Digital Video Endoscopy System offering classic features supporting the user for precission diagnosis & detection.

Key Features



HD Imaging



PIP/POP



Extreme Close up Endoscopy (eCUE)



Multiple O/P & I/P Video Signal



Compact & Light Weight

Specifications:

Power Supply	Voltage	100-240 VAC
Size	Voltage fluctuation	Within+10%
5,20	Frequency	50/60 Hz
	Frequency fluctuation	Within+3 Hz
	Consumption electric power	72VA
	Fuse rating	4 A, 250 V
	Fuse size	ø5.0 x 20.0 mm
	Dimensions & Weight	320(W) x 500(L) x 152(H) mm 10.3kg
Classification (medical	Type of protection against electric shock.	Class 1
electrical equipment)	Degree of protection against electric stock of applied part.	Type BF applied part. Where no classification mark appears, the device is a TYPE BF applied part.
	Degree of protection against explosion.	The Video Processor should be kept away from fammable gases.
Observation	Video signal output	VBS composite (PAL), Y/C, HDMI, HD-SDI, DVI-D Simultaneous output is possible.
	Video signal input	VBS composite (PAL), HDMI
	White balance	By pressing the white balance switch on the front panel, automatic white balance can be performed.
	Color tone adjustment	The following color tone adjustments may be using the color tone adjustment on the "System Setup" menu * "R" control:±10 steps * "B" control:±10 steps.
	Iris area selection	The area over which light is measured can be changed according the body part that is beir observed. The following to choices can be made using the iris mode selector on the front pareness. Normal observation Peak : When focusing on and/or observing a small bright a
	Image Enhancement mode setting	Structure enhancement in the endoscopic images can be enhanced electrically to increase the image sharpness.
	Gamma Correction	Gamma correction adjustments can be performed by adjusting the gamma value.
	Switching the enhancement modes	The enhancement level can be selected from 3 levels enhancement modes (Low, Med. and High) using theimage enhancement mode on the front button panel.
	Freeze screen display	The endoscopic image can be displayed as a frozen image, using the freeze switch key from endoscope.
	Brightness Selector	The brightness can be adjusted using the brightness adjustment on the "system setup/camera adjustment" menu.
	Remote control	The following function can be controlled by the endoscope's remote switches (Freeze, Recu Light, White Balance, mBLU, Enhancement, Contrast, Focus, Iris, Image Size, Remove Dat EXT In)
	mBLU	This mode is used to switch the normal mode to mBLU mode with white or blue light or vice versa.
De au ma ambabi	External Storage	USB Memory Stick
Documentation	Patient Data	The following data and modes can be displayed on the video monitor using the keyboard-ID number, patient Name, Sex, Age, Date of Birth, & Comments.
	Display the record state	Date/time (built in clock)
	Pre-procedural patient data	For a maximum of 40 patients- ID number, patient Name, Sex, Age,& Date of Birth
Image Storage & Retrieval	Memorization of selected setting	The following setting on the front panel are stored even when the video processor is turne OFF- System Setup, Setting the image record, Patient Data, Colour Tone,& White Balance
Memory Backup	Lithium Battery	Life: 5 years
Illumination	LED (White)	At Endoscope Tip
Air Feeding	Pump	Diaphragm type pump
Method	Airflow	"Low", "High", "Stop"
	Airflow Pressurization of water container (detachable)	Air pressurization of water container (detachable)

19 20 Index

^{*}Compatible with all 2H Series Endoscopes





HD VIDEO GASTROSCOPE

With multiple beneficial LED At Tip Technology



Model: OVG-2H (Compatible Video Processor: OVP-SEII Plus, OVP-SEII)

21

Key Features



Fiberless LED-At-Tip



Auxiliary Water Jet Channel




HD Imaging



Single Step, Waterproof Connector









Specifications:


Optical System	Field Of View	140°
optical system	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
	Optical Technology - mBLU	mBLU Mode 1
Insertion Tube	Distal End Outer Diameter	9.8 mm
	Insertion Tube Outer Diameter	9.5 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W)
		3. Objective Lens
		4. Instrument Channel Outlet
		5. Auxiliary Water Channel
Instrument Channel	Working Length	1055 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		30 cm 3/s
Bending Section	Angulation Range	UP 210°
		DOWN 90°
		RIGHT 100°
		LEFT 100°
Total Length		1365 mm
Compatible Video Processor		OVP-SEII Plus, OVP-SEII





ULTRA SLIM HD VIDEO GASTROSCOPE

With multiple beneficial LED At Tip Technology



Model: OVG-2HUS (Compatible Video Processor: OVP-SEII Plus, OVP-SEII)

Key Features



Fiberless LED-At-Tip



HD Imaging



Single Step, Waterproof Connector

23









Specifications:

Optical System	Field Of View	140°
option system	Direction Of View	0° Forward Viewing
	Depth Of Field	3-100 mm
Insertion Tube	Distal End Outer Diameter	7.8 mm
	Insertion Tube Outer Diameter	7.8 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W)
		3. Objective Lens
		4. Instrument Channel Outlet
Instrument Channel	Working Length	1055 mm
	Channel Inner Diameter	2.2 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		30 cm 3/s
Bending Section	Angulation Range	UP 210°
		DOWN 90°
		RIGHT 100°
		LEFT 100°
Total Length		1365 mm
Compatible Video Processor		OVP-SEII Plus, OVP-SEII





HD VIDEO COLONOSCOPE

New Tactile Insertion Tube (TIT) technology



Model: OVC-2H (L/I/S) (Compatible Video Processor: OVP-SEII Plus, OVP-SEII)

Key Features







HD Imaging



Tactile Insertion Tube



Auxiliary Water Jet Channel



Single Step, Waterproof Connector









Optical System	Field Of View	140°
Optical System	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
	Optical Technology - mBLU	mBLU Mode 1
Insertion Tube	Distal End Outer Diameter	13 mm
	Insertion Tube Outer Diameter	12.8 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W/B)
		3. Objective Lens
		4. Instrument Channel Outlet
		5. Auxiliary Water Channel
	Working Length	L: 1680 mm
		I: 1330 mm
		S: 730 mm
Instrument Channel	Channel Inner Diameter	3.7 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		25 cm %
Bending Section	Angulation Range	UP 180°
		DOWN 180°
		RIGHT 160°
		LEFT 160°
Total Length		L: 1980 mm
		I: 1630 mm
		S: 1030 mm
Compatible Video Processor		OVP-SEII Plus, OVP-SEII





HD VIDEO COLONOSCOPE

HD Paediatric Video Endoscope



Model: OVC-P2H (L/I) (Compatible Video Processor: OVP-SEII Plus, OVP-SEII)

Key Features







HD Imaging

27



Tactile Insertion Tube



Auxiliary Water Jet Channel



Single Step, Waterproof Connector









Specifications:

Optical System	Field Of View	140°	
Optical System	Direction Of View	Forward Viewing	
	Depth Of Field	3-100 mm	
	Optical Technology - mBLU	mBLU Mode 1	
Insertion Tube	Distal End Outer Diameter	11.8 mm	
	Insertion Tube Outer Diameter	11.8 mm	
	Distal End Enlarged	1. Air/Water Nozzle	
		2. Illumination LEDs (W/B)	
		3. Objective Lens	
		4. Instrument Channel Outlet	
		5. Auxiliary Water Channel	
	Working Length	L: 1680 mm	
		I: 1330 mm	
Instrument Channel	Channel Inner Diameter	3.2 mm	
	Minimum Visible Distance	3 mm from the distal end	
	Direction from which endo-therapy accessories enter and exit the endoscopic image.		
Air Flow Rate		25 cm 3/s	
Bending Section	Angulation Range	UP 180°	
		DOWN 180°	
		RIGHT 160°	
		LEFT 160°	
Total Length		L: 1980 mm	
		I: 1630 mm	
Compatible Video Processor		OVP-SEII Plus, OVP-SEII	





HD VIDEO DUODENOSCOPE

Excellent image quality & easy operability in Therapeutic procedure



Model: OVD-2H (Compatible Video Processor: OVP-SEII Plus, OVP-SEII)

29

Key Features



Fiberless LED-At-Tip



Single Step, Waterproof Connector









Specifications:

Optical System	Field Of View	100°	
Optical System	Direction Of View	Backward side viewing 5°	
	Depth Of Field	5-60 mm	
	Optical Technology - mBLU	mBLU Mode 1	
Insertion Tube	Distal End Outer Diameter	13.5 mm	
	Insertion Tube Outer Diameter	12.5 mm	
	Distal End Enlarged 3	1. Air/Water Nozzle	
	2 2 4	2. Objective Lens	
	15	3. Illumination LEDs (W)	
		4. Elevator	
		5. Instrument Channel Outlet	
Instrument Channel	Working Length	1260 mm	
	Channel Inner Diameter	4.2 mm	
	Minimum Visible Distance	10 mm from the Distal End	
	Direction from which endo-therapy accessories enter and exit the endoscopic image.		
Air Flow Rate		25 cm 3/s	
Bending Section	Angulation Range	UP 120°	
		DOWN 90°	
		RIGHT 110°	
		LEFT 90°	
Total Length		1545 mm	
Compatible Video Processor		OVP-SEII Plus, OVP-SEII	





Index



SmartEye-S3

Model: OVP-S3

Multi-featured Video Endoscopy System

Key Features



Fiberless LED-At-Tip



Multiple Video Signals





8 IN 1 Scope Connectivity



Scope Identification



Video Recording



Single Step, Waterproof Connector



Image Capturing



Compact & Light Weight

Specifications:

	Voltage	100-240 VAC
ower Supply	Voltage	
	Voltage fluctuation	Within ±10%
	Frequency	50/60 Hz
	Frequency fluctuation	Within ±3 Hz
	Consumption electric power	40 VA
	Fuse rating	4 A, 250 V
	Fuse size	ø5.0 x 20.0 mm
Size	Dimensions	445 (W) x 320 (L) x 152 (H) mm
	Weight	8.9 kg
	Type of protection against electric shock.	Class 1
Classification	Degree of protection against electric	Depends on the applied part. (The degree of protection
medical electrical	shock of applied part.	against electric shock of this product is BF type if the
equipment)	Shock of applied part.	mounting part be connected to this product is BF type. No
		that CF type is not compatible with this product).
	Degree of protection against explosion.	The Video Processor should be kept away
		from flammable gases.
Observation	Analog SD TV signal output	$Y/C \times 1$, VBS $\times 1$ simultaneous output is possible.
	Digital signal output	HDMI x 2 , DVI-D x 2 simultaneous output is possible.
	White balance	By pressing the white balance switch on the front
	Write Dalatice	
		panel, automatic white balance can be performed.
	Color tone adjustment	The following color tone adjustments may be using
		the color tone adjustment on the "System Setup"
		menu.* "R" control: ±5 steps, * "G" control: ±5
		steps,* "B" control: ±5 steps, * Chroma : ±5 steps.
	Imaga cattings	Contrast: ±5 steps, Brightness: ±5 steps, Gamma: ±5 st
	Image settings	
	Image size option	Small, Medium, Default & Large
	Light target	ALC: ±5 steps
	Iris area selection	The area over which light is measured can be changed
		according the body part that is being observed. The follow
		two choices can be made using the iris mode selector on
		keyboard and endoscope remote switch. Average: Norma
		·
		observation.Peak: When focusing on and/or observing a
		small bright area.
	Image enhancement mode setting	Fine patterns or edges in the endoscopic images can be
		enhanced electrically to increase the image sharpness.
	Gamma correction	Gamma correction adjustments can be performed by
	24	adjusting the gamma value.
	Conitabina the contractor	
	Switching the enhancement	The enhancement level can be selected from 3 Level
	modes	enhancement mode (Low, Med and High) using the image
		enhancement mode on the front button panel.
	Brightness selector	The brightness can be adjusted using the brightness
	3	adjustment on the "System & Image setting menu.
	Contract coloctor	
	Contrast selector	The contrast can be adjusted using the contrast adjustment of the contrast
		on the "System Image setting" menu.
	Zoom	1.3x
	Remote control	The following functions can be controlled by the endosco
		remote switches (Light, White Balance, Freeze, Zoom, Re
		Video, Iris, Record Image, Enhancement, Image Size, ml
		Remove Data, No Option)
	*mBLU (Optional feature)	This mode is used to switch the normal mode to mBLU
		mode - 1 & mBLU mode - 2
	Freeze screen display	The endoscopic image can be displayed as a frozen image
	1 reeze sereem dispilay	using the freeze switch key from endoscope.
		, ,
	Pre-Freeze	While capturing a image, the function analyzes previous
		images & display the sharpest image from the view.
Portable Memory	Media	USB mass storage
,	Image recording	PNG Format
	Movie recording	AVI
Documentation	Patient data	The following data and modes can be displayed on the vi
		monitor using the keyboard - ID number, Doctor Name,
		Patient Name, Sex, Age, Date of Birth & Comments.
	Display the record state	Date/time (built-in clock)
	Scope information data	The following scope information can be displayed on the
	Scope information data	
		monitor - Scope Model No., Serial No., Insertion Tube
		Diameter, Channel Diameter, Distal End Diameter &
		User details
USB Ports	Total 3 nos.	Keyboard & External Storage
Memory backup	Number of user settings	User settings up to 5nos. users
,	User function settings	Settings are held in video processor memory, even after
		video system is turned off.
	Lithium battery	
Till	·	Life: 5 years
Illumination	Fibreless LED (White)	At the endoscope's Tip
	Pump	Diaphragm type pump
	Airflow	3-mode (OFF, Low, High)
		Air pressurization of water container (detachable)
Water feeding		

*Compatible with all S2+ Series Endoscopes 31 32





VIDEO GASTROSCOPE

S2+ Series



Model: OVG-S2+ (Compatible Video Processor: OVP-S3)

33

Key Features



Fiberless LED-At-Tip



Remote Button



CMOS Sensor



Single Step, Waterproof Connector









Specifications:

Optical System	Field Of View	120°
option system	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
Insertion Tube	Distal End Outer Diameter 9.8 mm	
	Insertion Tube Outer Diameter	9.5 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W)
		3. Objective Lens
		4. Instrument Channel Outlet
Instrument Channel	Working Length	1055 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		30 cm 3/s
Bending Section	Angulation Range	UP 210°
		DOWN 90°
		RIGHT 100°
		LEFT 100°
Total Length		1365 mm
Compatible Video Processor		OVP-S3





ULTRA SLIM VIDEO GASTROSCOPE

S2+ Series



Model: OVG-S2+US (Compatible Video Processor: OVP-S3)

35

Key Features



Fiberless LED-At-Tip



Remote Button



CMOS Sensor











Specifications:


Optical System	Field Of View	120°
optical system	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
Insertion Tube	Distal End Outer Diameter	7.8 mm
	Insertion Tube Outer Diameter	7.8 mm
	Distal End Enlarged	Air/Water Nozzle Left Right Objective Lens Down Instrument Channel Outl
Instrument Channel	Working Length	1055 mm
	Channel Inner Diameter	2.2 mm
	Minimum Visible Distance	3.5 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 210°
		DOWN 90°
		RIGHT 100°
		LEFT 100°
Total Length		1365 mm
Compatible Video Processor		OVP-S3





VIDEO COLONOSCOPE

S2+ Series



Model: OVC-S2+
(Compatible Video Processor: OVP-S3)

37

Key Features



Fiberless LED-At-Tip




Remote Button



CMOS Sensor



Single Step, Waterproof Connector










Optical System	Field Of View	120°
	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
Insertion Tube	Distal End Outer Diameter	13 mm
	Insertion Tube Outer Diameter	12.8 mm
	Distal End Enlarged	1. Air/Water Nozzle
		2. Illumination LEDs (W)
		3. Objective Lens
		4. Instrument Channel Outlet
Instrument Channel	Working Length	1650 mm
	Channel Inner Diameter	3.7 mm
	Minimum Visible Distance	5 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Air Flow Rate		25 cm 3/s
Bending Section	Angulation Range	UP 180°
		DOWN 180°
		RIGHT 160°
		LEFT 160°
Total Length		1980 mm
Compatible Video Processor		OVP-S3





VIDEO COLONOSCOPE

HD Paediatric Video Endoscope S2+ Series



Model: OVC-PS2+L (Compatible Video Processor: OVP-S3)

39

Key Features



Fiberless LED-At-Tip




Remote Button



CMOS Sensor



Single Step, Waterproof Connector









Specifications:

Optical System	Field Of View	120°	
optical system	Direction Of View	Forward Viewing	
	Depth Of Field	3-100 mm	
Insertion Tube	Distal End Outer Diameter	11.8 mm	
	Insertion Tube Outer Diameter	11.8 mm	
	Distal End Enlarged	1. Air/Water Nozzle	
		2. Illumination LEDs (W)	
		3. Objective Lens	
		4. Instrument Channel Outlet	
Instrument Channel	Working Length	1680 mm	
	Channel Inner Diameter	3.2 mm	
	Minimum Visible Distance	3 mm from the Distal End	
	Direction from which endo-therapy accessories enter and exit the endoscopic image.		
Air Flow Rate		25 cm 3/s	
Bending Section	Angulation Range	UP 180°	
		DOWN 180°	
		RIGHT 160°	
		LEFT 160°	
Total Length		1980 mm	
Compatible Video Processor		OVP-S3	





VIDEO BRONCHOSCOPE

S2+ Series



Model: OVB-S2+ / OVB-XS2+ / OVB-PS2+ / OVB-XPS2+ / OVB-CS2+ (Compatible Video Processor: OVP-S3)

41

Key Features



Fiberless LED-At-Tip




Laser Compatible



CMOS Sensor



Single Step, Waterproof Connector









Specifications:

MODEL		OVB-S2+	OVB-XS2+	OVB-PS2+	OVB-XPS2+	OVB-CS2+
Optical System	Field of view	120°	120°	120°	120°	120°
	Direction of view	Forward View	Forward View	Forward View	Forward View	Forward View
	Depth of field	3-60mm	3-60mm	3-60mm	3-60mm	3-60mm
Insertion Tube	Distal End Outer Diameter	6 mm	6 mm	5 mm	5 mm	3.4 mm
	Insertion Tube Outer Diameter	6 mm	6 mm	5 mm	5 mm	3.3 mm
	Distal End Enlarged	1. Illumination				
		LEDs (W)				
		2. Image Sensor	2. Image Sensor	2. Image Sensor	2. Image Sensor	2. Image Sens
		3. Instrument				
		Channel	Channel	Channel	Channel	Channel
		Outlet	Outlet	Outlet	Outlet	Outlet
	Working Length	600 mm				
Instrument Channel	Channel Inner diameter	2.8 mm	3.2 mm	2 mm	2.2 mm	1.2 mm
	Minimum Visible Distance	3 mm from the				
		Distal End				
	Direction from which endo-					
	therapy accessories enter					
	and exit the endoscopic					
	image					
Bending Section	Angulation Range	UP 180°,				
		DOWN 130°				
Total Length		925 mm	925 mm	925 mm	925 mm	955 mm
Laser Compatibility		Nd: YAG,				
		810nm Diode				
Compatible Video Process	sor	OVP-S3	OVP-S3	OVP-S3	OVP-S3	OVP-S3





VIDEO PLEURASCOPE

A Slim Size Thoracoscope with Fiberless LED at Tip technology S2+ Series



Model: OTS-S2+ (Compatible Video Processor: OVP-S3)

Key Features



Fiberless LED-At-Tip



Laser Compatible



CMOS Sensor



Single Step, Waterproof Connector










Optical System	Field Of View	120°
	Direction Of View	Forward Viewing
	Depth Of Field	3-60 mm
insertion Tube	Distal End Outer Diameter 6 mm	
	Rigid Section Outer Diameter	6 mm
	Maximum Insertion Section Outer Diameter	6.7 mm
	Distal End Enlarged	1. Illumination LEDs (W)
		2. Image Sensor
Instrument Channel		3. Instrument Channel Outlet
	Working Length	300 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 160°
		DOWN 130°
Total Length		595 mm
aser Compatiblity		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-S3







Distal End Tip

Light Weight Design

Flexible Angulation





VIDEO LARYNGOSCOPE

S2+ Series



Model: OVL-S2+ (Compatible Video Processor: OVP-S3)

Key Features



Fiberless LED-At-Tip



Laser Compatible




CMOS Sensor



Single Step, Waterproof Connector









Optical System	Field Of View	120°
	Direction Of View	Forward Viewing
	Depth Of Field	3-60 mm
nsertion Tube	Distal End Outer Diameter	4.9 mm
	Insertion Tube Outer Diameter	4.8 mm
	Distal End Enlarged	1. Illumination LEDs (W)
		2. Image Sensor
		3. Instrument Channel Outlet
Instrument Channel	Working Length	350 mm
	Channel Inner Diameter	2 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 130°
		DOWN 130°
Total Length		645 mm
Laser Compatiblity		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-S3





VIDEO URETEROSCOPE

S2+ Series



Model: OUR-S2 / OUR-S2+ (Compatible Video Processor: OVP-S3)

47

Key Features



Fiberless LED-At-Tip



Laser Compatible



CMOS Sensor



Single Step, Waterproof Connector









Specifications:

MODEL:		OUR-S2	OUR-S2+
Optical System	Field Of View	120°	120°
, , , , , , , , , , , , , , , , , , , ,	Direction Of View	Forward Viewing	Forward Viewing
	Depth Of Field	3-60 mm	3-60 mm
Insertion Tube	Distal End Outer Diameter	3.3 mm	2.85 mm Bullet Shape
	Insertion Tube Outer Diameter	3.3 mm	2.85 mm
	Distal End Enlarged	1. Illumination LEDs (W)	1. Illumination LEDs (W)
		2. Image Sensor	2. Image Sensor
		3. Instrument Channel Outlet	3. Instrument Channel Out
Instrument Channel	Working Length	670mm	670 mm
	Channel Inner Diameter	1.2 mm	1.2 mm
	Minimum Visible Distance	3 mm from the Distal End	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.		
Bending Section	Angulation Range	UP 270°	UP 270°
		DOWN 270°	DOWN 270°
Total Length		965 mm	965 mm
Laser Compatiblity		Nd: YAG, 810nm Diode	Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-S3	OVP-S3





VIDEO CYSTOSCOPE

S2+ Series



Model: OCY-S2+ (Compatible Video Processor: OVP-S3)

49

Key Features



Fiberless LED-At-Tip



Laser Compatible



CMOS Sensor



Single Step, Waterproof Connector









Optical System	Field Of View	120°	
	Direction Of View	Forward Viewing	
	Depth Of Field	3-60 mm	
Insertion Tube	Distal End Outer Diameter	2.80 mm Bullet Shape	
	Insertion Tube Outer Diameter	5 mm	
	Distal End Enlarged	1. Illumination LEDs (W)	
		2. Image Sensor	
		3. Instrument Channel Outlet	
Instrument Channel	Working Length	380 mm	
	Channel Inner Diameter	2 mm	
	Minimum Visible Distance	3 mm from the Distal End	
	Direction from which endo-therapy accessories enter and exit the endoscopic image.		
Bending Section	Angulation Range	UP 210°	
		DOWN 130°	
Total Length		675 mm	
Laser Compatiblity		Nd: YAG, 810nm Diode	
Compatible Video Processor		OVP-S3	











Insight-II

Model: OVP-IN2

Next Generation Digital Video Endoscopy Processor

COMPATIBLE WITH

• Bronchoscope • Pleurascope • Laryngoscope • Ureteroscope • Cystoscope • Cholangioscope

Key Features



Fiberless LED-At-Tip



Multiple Video Signals



6 IN 1 Scope Connectivity



Scope Identification



Video Recording



Image Capturing



Compact & Light Weight

*Compatible with all IN+ Series Endoscopes

Specifications:

Dower Supply	Voltage	100-240 VAC
Power Supply	Voltage fluctuation	Within ±10%
	Frequency	50/60 Hz
	Frequency fluctuation	Within ±3 Hz
	Consumption electric power	18 VA
	Fuse rating	3 A, 250 V
	Fuse size	ø5.0 x 20.0 mm
Size	Dimensions	300 (W) x 295 (L) x 80 (H) mm
	Weight	5 kg
	Type of protection against electric shock.	Class 1
Classification	Degree of protection against electric	Type BF applied part. Where no classification
(medical electrical	shock of applied part.	mark appears, the device is a type BF applied part.
equipment)	Degree of protection against explosion.	The Video Processor should be kept away
		from flammable gases.
Observation	Analog SD TV signal output	Y/C x 1, VBS x 1 simultaneous output is possible.
	Digital signal output	HDMI x 2 , DVI-D x 2 simultaneous output is possible.
	White balance	By pressing the white balance switch on the front
	Write balance	panel, automatic white balance can be performed.
	Color tone adjustment	The following color tone adjustments may be using
		the color tone adjustment on the "System Setup"
		menu.* "R" control: ±5 steps, * "G" control: ±5
		steps,* "B" control: ±5 steps, * Chroma : ±5 steps.
	Image settings	Contrast: ±5 steps, Brightness: ±5 steps, Gamma: ±5 steps
	Image settings Image size option	Small, Medium, Default & Large
	Light target	ALC: ±5 steps
	Iris area selection	The area over which light is measured can be changed
		according the body part that is being observed. The followi
		two choices can be made using the iris mode selector on
		keyboard and endoscope remote switch. Average: Normal
		observation.Peak: When focusing on and/or observing a
		, , , , , , , , , , , , , , , , , , , ,
		small bright area.
	Image enhancement mode setting	Fine patterns or edges in the endoscopic images can be
		enhanced electrically to increase the image sharpness.
	Gamma correction	Gamma correction adjustments can be performed by
		adjusting the gamma value.
	Switching the enhancement	The enhancement level can be selected from 3 Level
	modes	
	modes	enhancement mode (Low, Med and High) using the image
		enhancement mode on the front button panel.
	Brightness selector	The brightness can be adjusted using the brightness
		adjustment on the "System & Image setting menu.
	Contrast selector	The contrast can be adjusted using the contrast adjustmen
	Contrast selector	on the "System Image setting" menu.
	7	, 3 3
	Zoom	1.3x
	Remote control	Multiple programmable functions can be controlled by
		the endoscope's 2Nos.remote switches
	*mBLU (Optional feature)	This mode is used to switch the normal mode to mBLU
	indeo (optional reduite)	mode - 1 & mBLU mode - 2
	Fuere emerical districtions	
	Freeze screen display	The endoscopic image can be displayed as a frozen image,
		using the freeze switch key from endoscope.
	Pre-Freeze	While capturing a image, the function analyzes previous
		images & display the sharpest image from the view.
Portable Memory	Media	USB mass storage
•	112212	
Documentation	Patient data	The following data and modes can be displayed on the vide
		monitor using the keyboard * ID number * Doctor Name *
		Patient Name * Sex, Age * Date of Birth * Comments.
	B: 1 11 1 1 1 1	
	Display the record state	* Date/time (built-in clock)
	Scope information data	The following scope information can be displayed on the
		monitor * Scope Model No. * Serial No. * Insertion Tube
		Diameter * Channel Diameter * Distal End Diameter
		* User details
USB Ports	Total 3 nos.	Keyboard & External Storage
Memory backup	Number of user settings	User settings up to 5nos. users
-	User function settings	Settings are held in video processor memory, even after
		video system is turned off.
	Lithium battery	Life: 5 years
Illumination	Fibreless LED (White)	At the endoscope's Tip





VIDEO BRONCHOSCOPE

Insight Plus Series



Model: OVB-IN+ / OVB-XIN+ / OVB-PIN+ / OVB-XPIN+ / OVB-CIN+ (Compatible Video Processor: OVP-IN2)

53

Key Features



Fiberless LED-At-Tip



Laser Compatible




CMOS Sensor











Specifications:

MODEL: OCY-	S2+	OVB-IN+	OVB-XIN+	OVB-PIN+	OVB-XPIN+	OVB-CIN+
Optical System	Field Of View	120°	120°	120°	120°	120°
. ,	Direction Of View	Forward Viewing				
	Depth Of Field	3-60 mm				
Insertion Tube	Distal End Outer Diameter	6 mm	6 mm	5 mm	5 mm	3.4 mm
	Insertion Tube Outer Diameter	6 mm	6 mm	5 mm	5 mm	3.3 mm
	Distal End Enlarged	1. Illumination LEDs (W)	1. Illumination LEDs (W)	1. Illumination LEDs (W)	1. Illumination LEDs (W)	1. Illumination LEDs (W)
		2. Image Sensor	2. Image Sensor	2. Image Sensor	2. Image Sensor	2. Image Senso
		3. Instrument Channel Outlet				
	Working Length	600 mm				
Instrument Channel	Channel Inner Diameter	2.8 mm	3.2 mm	2 mm	2.2 mm	1.2 mm
	Minimum Visible Distance	3 mm from the Distal End				
	Direction from which endo-therapy accessories enter and exit the endoscopic image.					
Bending Section	Angulation Range	UP 180°				
		DOWN 130°				
Total Length		925 mm	925 mm	925 mm	925 mm	955 mm
Laser Compatiblity		Nd: YAG, 810nm Diode				
Compatible Video Processor		OVP-IN2	OVP-IN2	OVP-IN2	OVP-IN2	OVP-IN2





VIDEO PLEURASCOPE

Insight Plus Series
A Slim Size Thoracoscope with Fiberless LED at Tip technology



Model: OTS-IN+
(Compatible Video Processor: OVP-IN2)

55

Key Features



Fiberless LED-At-Tip



Laser Compatible



CMOS Sensor



Single Click Connector










Specifications:

ptical System	Field Of View	120°
-	Direction Of View	Forward Viewing
	Depth Of Field	3-60 mm
Insertion Tube	Distal End Outer Diameter	6 mm
	Rigid Section Outer Diameter	6 mm
	Maximum Insertion Section Outer Diameter	6.7 mm
	Distal End Enlarged	1. Illumination LEDs (W)
		2. Image Sensor
Instrument Channel		3. Instrument Channel Outlet
	Working Length	300 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 160°
		DOWN 130°
Total Length		595 mm
Laser Compatiblity		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-IN2





Distal End Tip

Light Weight Design

Flexible Angulation





VIDEO LARYNGOSCOPE

Insight Plus Series



Model: OVL-IN+
(Compatible Video Processor: OVP-IN2)

57

Key Features



Fiberless LED-At-Tip



Laser Compatible




CMOS Sensor



Single Click Connector









Optical System	Field Of View	120°
	Direction Of View	Forward Viewing
	Depth Of Field	3-60 mm
Insertion Tube	Distal End Outer Diameter	4.9 mm
	Insertion Tube Outer Diameter	4.8 mm
	Distal End Enlarged	1. Illumination LEDs (W)
		2. Image Sensor
		3. Instrument Channel Outlet
Instrument Channel	Working Length	350 mm
	Channel Inner Diameter	2 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 130°
		DOWN 130°
Total Length		645 mm
Laser Compatiblity		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-IN2





VIDEO URETEROSCOPE

Insight Plus Series



Model:OUR-IN+
(Compatible Video Processor: OVP-IN2)

59


Key Features



Fiberless LED-At-Tip



Laser Compatible



CMOS Sensor









MODEL:		OUR-IN+
Optical System	Field Of View	120°
7.55	Direction Of View	Forward Viewing
	Depth Of Field	3-60 mm
Insertion Tube	Distal End Outer Diameter	2.85 mm Bullet Shape
	Insertion Tube Outer Diameter	2.85 mm
	Distal End Enlarged	1. Illumination LEDs (W)
		2. Image Sensor
		3. Instrument Channel Outlet
Instrument Channel	Working Length	670 mm
	Channel Inner Diameter	1.2 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 270°
		DOWN 270°
Total Length		965 mm
Laser Compatiblity		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-IN2





VIDEO CYSTOSCOPE

Insight Plus Series



Model: OCY-IN+
(Compatible Video Processor: OVP-IN2)

61

Key Features



Fiberless LED-At-Tip



Laser Compatible



CMOS Sensor




Single Click Connector









Optical System	Field Of View	120°
, , , , , , , , , , , , , , , , , , , ,	Direction Of View	Forward Viewing
	Depth Of Field	3-60 mm
Insertion Tube	Distal End Outer Diameter	2.80 mm Bullet Shape
	Insertion Tube Outer Diameter	5 mm
	Distal End Enlarged	1. Illumination LEDs (W)
		2. Image Sensor
		3. Instrument Channel Outlet
Instrument Channel	Working Length	380 mm
	Channel Inner Diameter	2 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 210°
		DOWN 130°
Total Length		675 mm
Laser Compatiblity		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-IN2





VIDEO CHOLANGIOSCOPE

Ultra slim baby endoscope Insight Plus Series



Model: OBS-IN+ (Compatible Video Processor: OVP-IN2)

Key Features



Fiberless LED-At-Tip



CMOS Sensor





of view






120° wide field Laser Compatible


Single Click Connector







Optical System	Field Of View	120°
Optical System	Direction Of View	Forward Viewing
	Depth Of Field	2-60 mm
Insertion Tube	Distal End Outer Diameter	3.3 mm
	Insertion Tube Outer Diameter	3.4 mm
	Distal End Enlarged	1. Illumination LEDs (W)
		2. Image Sensor
		3. Instrument Channel Outlet
Instrument Channel	Working Length	1920 mm
	Channel Inner Diameter	1.2 mm
	Minimum Visible Distance	2 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 75°
		DOWN 90°
Total Length		2215 mm
Laser Compatiblity		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-IN2







Distal End Tip

Light Weight Design

Flexible Angulation













SmartEye-S3

Model: OVP-S3

Multi-featured Video Endoscopy System

Key Features



Fiberless LED-At-Tip



Multiple Video Signals





Connectivity



Scope Identification



Video Recording



Image Capturing







Compact & Light Weight

Specifications:

	Voltage	100-240 VAC
Power Supply		
	Voltage fluctuation	Within ±10%
	Frequency	50/60 Hz
	Frequency fluctuation	Within ±3 Hz
	Consumption electric power	40 VA
	Fuse rating	4 A, 250 V
	Fuse size	ø5.0 x 20.0 mm
Size	Dimensions	445 (W) x 320 (L) x 152 (H) mm
JIZC .	Weight	8.9 kg
	Type of protection against electric shock.	Class 1
Classification	Degree of protection against electric	Depends on the applied part. (The degree of protection
medical electrical	shock of applied part.	
equipment)	Shock of applied part.	against electric shock of this product is BF type if the
equipment)		mounting part be connected to this product is BF type. Note
		that CF type is not compatible with this product).
	Degree of protection against explosion.	The Video Processor should be kept away
		from flammable gases.
Observation	Analog SD TV signal output	Y/C x 1, VBS x 1 simultaneous output is possible.
22301 1441011	Digital signal output	HDMI x 2 , DVI-D x 2 simultaneous output is possible.
	White balance	
	white balance	By pressing the white balance switch on the front
		panel, automatic white balance can be performed.
	Color tone adjustment	The following color tone adjustments may be using
		the color tone adjustment on the "System Setup"
		menu.* "R" control: ±5 steps, * "G" control: ±5
		steps,* "B" control: ±5 steps, * Chroma : ±5 steps.
	Image cettings	
	Image settings	Contrast: ±5 steps, Brightness: ±5 steps, Gamma: ±5 step
	Image size option	Small, Medium, Default & Large
	Light target	ALC: ±5 steps
	Iris area selection	The area over which light is measured can be changed
		according the body part that is being observed. The following
		two choices can be made using the iris mode selector on
		keyboard and endoscope remote switch. Average: Normal
		observation.Peak: When focusing on and/or observing a
		small bright area.
	Image enhancement mode setting	Fine patterns or edges in the endoscopic images can be
		enhanced electrically to increase the image sharpness.
	Gamma correction	Gamma correction adjustments can be performed by
	Carrina correction	adjusting the gamma value.
	0 11 11 1	
	Switching the enhancement	The enhancement level can be selected from 3 Level
	modes	enhancement mode (Low, Med and High) using the image
		enhancement mode on the front button panel.
	Brightness selector	The brightness can be adjusted using the brightness
	J	adjustment on the "System & Image setting menu.
	Combined collection	
	Contrast selector	The contrast can be adjusted using the contrast adjustment
		on the "System Image setting" menu.
	Zoom	1.3x
	Remote control	The following functions can be controlled by the endoscopes
	Remote control	remote switches (Light, White Balance, Freeze, Zoom, Reco
		, , ,
		Video, Iris, Record Image, Enhancement, Image Size, mBLU
		Remove Data, No Option)
	*mBLU (Optional feature)	This mode is used to switch the normal mode to mBLU
	moto (optional reactive)	mode - 1 & mBLU mode - 2
	Former common Prof.	
	Freeze screen display	The endoscopic image can be displayed as a frozen image,
		using the freeze switch key from endoscope.
	Pre-Freeze	While capturing a image, the function analyzes previous
		images & display the sharpest image from the view.
Portable Momen	Media	USB mass storage
Portable Memory		
	Image recording	PNG Format
	Movie recording	AVI
Documentation	Patient data	The following data and modes can be displayed on the video
_ commentation	1 dilette data	monitor using the keyboard * ID number * Doctor Name *
		Patient Name * Sex, Age * Date of Birth * Comments.
	Display the record state	* Date/time (built-in clock)
	Scope information data	The following scope information can be displayed on the
		monitor * Scope Model No. * Serial No. * Insertion Tube
		Diameter * Channel Diameter * Distal End Diameter
LICD D t.	Tabel 2 and	* User details
USB Ports	Total 3 nos.	Keyboard & External Storage
Memory backup	Number of user settings	User settings up to 5nos. users
	User function settings	Settings are held in video processor memory, even after
		video system is turned off.
	Lithium battery	Life: 5 years
Illumination	Fibreless LED (White)	At the endoscope's Tip
munnination		·
	Pump	Diaphragm type pump
	Airflow	3-mode (OFF, Low, High)
Water feeding	Airflow	3-mode (OFF, Low, High) Air pressurization of water container (detachable)


*Compatible with all S2+ Series Endoscopes





VIDEO VETERINARY ENDOSCOPE

S2+ Series



Model: OVV-S2+ (I/S) (Compatible Video Processor: OVP-S3)

69

Key Features



Fiberless LED-At-Tip




Remote Button



CMOS Sensor



Single Step, Waterproof Connector









Specifications:

Optical System	Field Of View	120°
optical system	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
Insertion Tube	Distal End Outer Diameter	9.8 mm
	Insertion Tube Outer Diameter	9.8 mm
	Distal End Enlarged	Air/Water Nozzle Up Illumination LEDs (V Right Objective Lens Down Instrument Channel Outl
Instrument Channel	Working Length	I-1400 mm/ S- 1055 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 180°
		DOWN 14°
		RIGHT 120°
		LEFT 120°
Total Length		I-1710 mm/ S-1355mm
Compatible Video Processor		OVP-S3





VIDEO VETERINARY ENDOSCOPE

S2+ Series



Model: OVV-S2+US (I/S) (Compatible Video Processor: OVP-S3)

71

Key Features



Fiberless LED-At-Tip




Remote Button



CMOS Sensor



Single Step, Waterproof Connector









Specifications:

Optical System	Field Of View	120°
optical system	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
Insertion Tube	Distal End Outer Diameter	7.8 mm
	Insertion Tube Outer Diameter	7.8 mm
	Distal End Enlarged	Air/Water Nozzle Up Illumination LEDs (W) Left Right Objective Lens Down Instrument Channel Outlet
Instrument Channel	Working Length	I-1400 mm/ S- 1055 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 180°
		DOWN 120°
		RIGHT 120°
		LEFT 120°
Total Length		I-1710 mm/ S-1355mm
Compatible Video Processor		OVP-S3







SmartEye-II

Model: OVP-SEII

High Definition Video Endoscopy System

SmartEye-II is a High Definition Digital Video Endoscopy System offering classic features supporting the user for precission diagnosis & detection.

Key Features



HD Imaging



PIP/POP



Simultaneous Video Output



Compact & Light Weight

Specifications:


Power Supply	Voltage	100-240 VAC
Size	Voltage fluctuation	Within+10%
	Frequency	50/60 Hz
	Frequency fluctuation	Within+3 Hz
	Consumption electric power	72VA
	Fuse rating	4 A, 250 V
	Fuse size	ø5.0 x 20.0 mm
	Dimensions & Weight	320(W) x 500(L) x 152(H) mm 10.3kg
Classification (medical	Type of protection against electric shock.	Class 1
electrical equipment)	Degree of protection against electric stock of applied part.	Type BF applied part. Where no classification mark appears, the device is a TYPE BF applied part.
	Degree of protection against explosion.	The Video Processor should be kept away from fammable gases.
Observation	Video signal output	VBS composite (PAL), Y/C, HDMI, HD-SDI, DVI-D Simultaneous output is possible.
	Video signal input	VBS composite (PAL), HDMI
	White balance	By pressing the white balance switch on the front panel, automatic white balance can be performed.
	Color tone adjustment	The following color tone adjustments may be using the color tone adjustment on the "System Setup" menu * "R" control:±10 steps * "B" control:±10 steps.
	Iris area selection	The area over which light is measured can be changed according the body part that is being observed. The following to choices can be made using the iris mode selector on the front pa Average: Normal observation Peak: When focusing on and/or observing a small bright are
	Image Enhancement mode setting	Structure enhancement in the endoscopic images can be enhanced electrically to increase the image sharpness.
	Gamma Correction	Gamma correction adjustments can be performed by adjusting the gamma value.
	Switching the enhancement modes	The enhancement level can be selected from 3 levels enhancement modes (Low, Med. and High) using theimage enhancement mode on the front button panel.
	Freeze screen display	The endoscopic image can be displayed as a frozen image, using the freeze switch key from endoscope.
	Brightness Selector	The brightness can be adjusted using the brightness adjustment on the "system setup/camera adjustment"menu.
	Remote control	The following function can be controlled by the endoscope's remote switches (Freeze, Recor Light, White Balance, mBLU, Enhancement, Contrast, Focus, Iris, Image Size, Remove Data EXT In)
	mBLU	This mode is used to switch the normal mode to mBLU mode with white or blue light or vice versa.
	External Storage	USB Memory Stick
Documentation	Patient Data	The following data and modes can be displayed on the video monitor using the keyboard-ID number,patient Name, Sex, Age, Date of Birth, Comments.
	Display the record state	Date/time (built in clock)
Image Storage	Pre-procedural patient data Memorization of selected setting	For a maximum of 40 patients- ID number, patient Name, Sex, Age, Date of Birth The following setting on the front panel are stored even when the video processor is turned
& Retrieval		OFF- System Setup, Setting the image record, Patient Data, Colour Tone, White Balance.
Memory Backup	Lithium Battery	Life: 5 years
Illumination	LED (White)	At Endoscope Tip
Air Feeding	Pump	Diaphragm type pump
Method	Airflow Airflow Pressurization of water container (detachable)	"Low", "High", "Stop" Air pressurization of water container (detachable)

^{*}Compatible with all 2H Series Endoscopes





VIDEO VETERINARY ENDOSCOPE



Model: OVV-2H+ (L/I/S) (Compatible Video Processor: OVP-SEII)

Key Features



Fiberless LED-At-Tip




HD Imaging



Single Step, Waterproof Connector

75









Specifications:

Optical System	Field Of View	140°
optical system	Direction Of View	Forward Viewing
	Depth Of Field	3-100 mm
Insertion Tube	Distal End Outer Diameter	9.8 mm
	Insertion Tube Outer Diameter	9.8 mm
	Distal End Enlarged	Air/Water Nozzle Up Illumination LEDs (W) Left Right Objective Lens Down Instrument Channel Outlet
Instrument Channel	Working Length	L-3000 mm/ I-1400 mm / S- 1055 mm
	Channel Inner Diameter	2.8 mm
	Minimum Visible Distance	3 mm from the Distal End
	Direction from which endo-therapy accessories enter and exit the endoscopic image.	
Bending Section	Angulation Range	UP 180°
		DOWN 140°
		RIGHT 120°
		LEFT 120°
Total Length		L-3310 mm/ I-1710 mm / S- 1355 mm
Compatible Video Processor		OVP-S3



















Model: OLS

Endoplus HD

Full HD Laproscopic Camera System

Clear, Crisp & Brighter Imaging



HD Imaging



CMOS

CMOS Sensor



Single Step, Water Proof Connector



Model: OVP-EP

Optical Zoom



Ottomed LED Light Source

High end Light Source providing natural colors of video of an endoscopy camera.

A long life LED lamp

Features

- Long life LED with (5000 Hrs)
- Brightness control level 50% to 100%.
- Variable temperature settings for heating protection or undesired heat generation
- Flicker free steady LED light source.

Specification

- Ambient Temperature 10°c 40°c
- Relative Humidity 30% 85%
- Input Voltage 100v 240v

79 80 Index













Endocarbon-IIIi


Endoscopy CO₂ Insufflator

The Intelligent, Efficient & Compact machine for Patient comfort.

Highly efficient Intelligent gas flow.

Advance micro-controller technology.

Special customized & pre-programmed single touch mode for POEM procedure.



83



Key Features

- · Actual Gas Input Pressure Indicator
- Light Weight
- Compact Size
- Universal compatibility
- Single touch flow selection
- Eco-friendly

Specifications:

	MODEL	- m.EC-IIIi
Operating Environment		
Ambient Temperature		10°c – 40°c
Relative Humidity		30% to 85%
Main Unit Specifications		
Pneumatics		
Gas		CO ₂ (Medical Grade)
Input Pressure		1 bar - 7 bar
Flow Rate	POEM	1.0 L/min ±10%
	LOW	3 L/min ±10%
	HIGH	4.5 L/min ±10%
Intelligent Flow Control		Applicable
Power Supply		
Input Voltage		100V - 240V
Allowable Voltage Fluctuation	1	±10%
Power Rating		24VA
Frequency		50Hz / 60Hz
Allowable Frequency Fluctuat	tion	±3 Hz
Dimensions		
Length		220mm
Width		170mm
Height		145mm
Weight		Approx. 2.3 Kg








Endocarbon-III

Endoscopy CO₂ Insufflator

Efficient & Compact machine for Patient comfort.

Advance micro-controller technology. Special customized & pre-programmed single touch mode for POEM procedure.



85

Gas Flow Rate Indicator Select Flow Button -Gas Start/Stop Button Easy to Use... CO₂ Gas Outlet

Key Features

- · Actual Gas Input Pressure Indicator
- Light Weight
- Compact Size
- Universal compatibilitySingle touch flow selection

Specifications:

	MODEL	- m.EC-III
Operating Environment		
Ambient Temperature		10°c - 40°c
Relative Humidity		30% to 85%
Main Unit Specifications		
Pneumatics		
Gas		CO ₂ (Medical Grade)
Input Pressure		1 bar - 7 bar
Flow Rate	POEM	1.0 L/min ±10%
	LOW	3.5 L/min ±10%
	HIGH	5.0 L/min ±10%
Power Supply Input Voltage		100V - 240V
Allowable Voltage Fluctuation		±10%
Power Rating		24VA
Frequency		50Hz / 60Hz
requericy		
	ion	±3 Hz
Allowable Frequency Fluctuat	ion	±3 Hz
Allowable Frequency Fluctuat Dimensions	ion	±3 Hz
Allowable Frequency Fluctuat Dimensions Length	ion	T
Allowable Frequency Fluctuat Dimensions Length Width Height	ion	155 mm







RECOMMENDED APPLICATIONS

- · Auxiliary water jet irrigation for washing colonic mucosa.
- Auxiliary water jet/Biopsy channel irrigation for gastric washing, bleed etc.

87

• Quick filling for Endoscopic Ultrasound procedure.



- Reverse flow technology for precision flow control.
- Universal Compatibility.
- In-built safety alarm.
- Led Indication.
- · Noise less operation.
- Easy selection of water flow rate.
- Foot switch control.

Specifications:

MOD	EL - m.EP-II
Operating Environment	
Ambient Temperature	10°c - 40°c
Relative Humidity	30% to 85%
Main Unit Specifications	,
Power Supply	
Input Voltage	100V - 240V
Allowable Voltage Fluctuation	±10%
Current	0.4A - 0.2A
Frequency	50Hz / 60Hz
Water Flow Range	100ml/min - 900ml/min
Foot Switch Type	Pedal Type
Flow Indicator	LED
Dimensions	
Length	255 mm
Width	265 mm
Height	150 mm
Weight	Approx. 2.7 Kg Without Water Tank







Key Features



Detergent wash to remove sticky residue on Endoscopes



Dual Leakage Testing Port



Always visible disinfection, alcohol, detergent solution level for easy top up!



Instant thermal paper print-outs of endoscope data with re-processing protocol. History can also be printed

Specifications:

Alcohol Flush

	MODEL -	OEW-3	
Ambient temperature		10°c – 40°c	
Relative Humidity		30% to 85%	
Water Supply Requirement	Flow Rate	12L/min – 16L/min	
	Pressure	0.2 MPa – 0.5 MPa	
Main Unit Specification			
No. of Accommodated Endos	scopes	2 nos.	
Disinfectant Tank Capacity		25L	
Detergent tank capacity		1.75 Litres	
Alcohol tank capacity		1.75 Litres	
Leakage Testing		Bubble detecting during immersion (Visual)	
Operating Modes			
Auto Mode			
Manual Mode	Operations	Time Range	
	Wash Time	1 – 10 min	
	Detergent wash time	1 – 99 min	
	Disinfectant time	1 – 99 min	
	Dry time	1 – 10 min	
Alarm		At Low level	
		End of process	
Display		Touch Screen Panel	
Dimension		750 x 690 x 1040 (L x W x H)	
Weight		80 Kg Approx	
Power Supply			
Input Voltage		220 V	
Allowable Voltage Fluctuation		±10%	
Input Current		1.5 A	
Frequency		50Hz	
Allowable Frequency Fluctua	tion	±3	

Qubey-III

Fully Automatic Twin Scope Endowasher System with Alcohol flushing feature.

Minimize your Hygiene & Safety concerns.

Disinfect your scope in the **Qubey-III Endowasher** after every procedure to maximize safety of the Patient, Staff and yourself.

89







360° Turbo Jet

Qubey-II

Greater Visiblity for improved Diagnosis & Detection

Fully Automatic Twin Scope Endowasher System 360° Auto Disc Rotation Water Jet

Minimize your Hygiene & Safety concerns.

Disinfect your scope in the **Qubey-II Endowasher** after every procedure to maximize safety of the Patient, Staff and yourself.

91

Key Features



Detergent wash to remove sticky residue on Endoscopes



Dual Leakage Testing Port



Always visible disinfection, detergent solution level for easy top up!



Instant thermal paper print-outs of endoscope data with re-processing protocol. History can also be printed



Touch Screen Interface
Simple menu based
feather-touch operation.
Spill-proof long life touch screen

Specifications:

	MODE	EL - OEW-2
Ambient temperature		10°c - 40°c
Relative Humidity		30% to 85%
Water Supply Requirement	Flow Rate	12L/min – 16L/min
water Supply Requirement	Pressure	0.2 MPa - 0.5 MPa
Main Unit Specification		
No. of Accommodated End	oscopes	2 nos.
Disinfectant Tank Capacity		17L
Detergent Tank Capacity		1.5L
Leakage Testing		Bubble detecting during immersion (Visual)
Operating Modes		
Auto Mode		
Manual Mode	Operations	Time Range
	Wash Time	1 - 10 min
	Disinfectant time	1 – 99 min
	Dry time	1 – 10 min
	Detergent wash time	1 - 10 min
Alarm		At Low level
		End of process
Display		Touch Screen Panel
Dimension		600 x 600 x 600mm (L x W x H)
Weight		60 Kg Approx
Power Supply		
Input Voltage		220 V
Allowable Voltage Fluctuation		±10%
Input Current		1.5 A
Frequency		50Hz
Allowable Frequency Fluctu	uation	±3







Qubey-T

Greater Visiblity for improved Diagnosis & Detection

Fully Automatic Twin Scope Endowasher System

Minimize your Hygiene & Safety concerns.

Disinfect your scope in the **Qubey-T Endowasher** after every procedure to maximize safety of the Patient, Staff and yourself.

93

Key Features



Detergent wash to remove sticky residue on Endoscopes



Leakage Testing Port



Always visible disinfection solution level for easy top up!



Instant thermal paper print-outs of endoscope data with re-processing protocol. History can also be printed



Touch Screen Interface
Simple menu based
feather-touch operation.
Spill-proof long life touch screen

Specifications:

Ambient temperature		10°c – 40°c	
Relative Humidity		30% to 85%	
Water Supply Requirement	Flow Rate	12L/min – 16L/min	
water Supply Requirement	Pressure	0.1 MPa - 0.4 MPa	
Main Unit Specification			
No. of Accommodated Endos	copes	2 nos.	
Disinfectant Tank Capacity		15L	
Leakage Testing		Bubble detecting during immersion (Visual)	
Operating Modes			
Auto Mode			
Manual Mode	Operations	Time Range	
	Wash Time	1 – 20 Cycles	
	Disinfectant time	1 – 99 min	
	Dry time	1 – 10 min	
Alarm		At Low level	
		End of process	
Display		Touch Screen Panel	
Memory Storage		50 Cycles	
Dimension		600 x 600 x 600mm (L x W x H)	
Weight		55 Kg Approx	
Power Supply			
Input Voltage		220 V	
Allowable Voltage Fluctuation Input Current		±10%	
		1.5 A	
Frequency		50Hz	
Allowable Frequency Fluctuation		±3	













End Hanger 2

For Endoscope Storage

Improper storage is the most common reason for endoscope damage. EndoHanger 2 is a safe & secure endoscope storage tool to safely and securely store your endoscopes. Proper storage of endoscopes also plays an important part in endoscope reprocessing to ensure Patient and Technician safety.

Key Features

- Wall mounted Endoscope Hanger: Endohanger 2 is designed to be easily mounted on the wall or in the endoscope cabinet.
- Safer & Secure storage: Endohanger 2 helps save on costly repair costs and downtime by allowing proper storage of your expensive endoscopes.

95

- **Material:** Made from long-lasting durable single piece Acrylic which is easy to clean and wipe.
- Universal Compatibility.
- Can be costumed designed for EUS Endoscopes.

EndOHanger

For Endoscope Storage

Improper storage is the most common reason for endoscope damage. EndoHanger 2 is a safe & secure endoscope storage tool to safely and securely store your endoscopes. Proper storage of endoscopes also plays an important part in endoscope reprocessing to ensure Patient and Technician safety.

Key Features

- · For easy & safe endoscope hanging.
- Universal compatibility.
- Special grade rust free strong material.
- · Smooth surfaces for scratch free hanging of endoscope.









Model: OVM-24HD-A

97

Full HD Medical Grade Monitor

With precision color reproduction and smooth display for Ottomed Endoscopy Systems



Rear Panel Side Panel

Top Panel



Bottom Panel

Key Features

1. Accurate Color Rendering

Mitra endoscopic operating room professional displays utilize 12BIT technology, which allows for over 68.7 billion colors to be shown linearly, thus accurately rendering the color of human tissue, enabling endoscopists to make more accurate judgements.

2. Excellent Layering

Mitra endoscopic operating room professional displays use medical LCD panels with the highest definition. The panels possess fine color, high brightness, and high contrast, thus clearly representing each layer of tissue in endoscopy in its proper order, providing clear images for detection and diagnosis.

3. Smooth and Dynamic Images

The response time of endoscopic operating room professional displays is less than 30ms, thus ensuring the smooth display of images from standard definition to high definition, as well as showing indications of hemorrhage and other important indicators during endoscopy, with remarkable clarity and without any ghosting or shaking.

4. Compliance with Operating Room Sterilization and Anti-interference Requirements

Mitra endoscopic operating room professional displays adopt a white, dust-free design, which complies with operating room

sterilization requirements for electronic equipment. The metal-encased design can effectively prevent dust and fluid contamination during surgery. The metal casing body and the built-in shielding circuit can effectively prevent interference to the monitor from electronic surgical equipment ensuring the stability of the display.

5. Multiple Interface Types

Mitra endoscopic operating room professional displays have the compatibility for multiple analog and digital signals, and are perfectly compatible and suitable with and fully meeting the demands of Ottomed Endoscopy systems.

Specifications:

MODEL : OVM-24HD-A		
Backlight	LED	
Size	23.8°	
Active Display	527.04 (H) x 296.46 (V) mm	
Resolution	1920x1080 (HD1080p)	
Pixel Pitch	0.270 (H) x 0.270 (V) mm	
Response Time	5ms	
Brightness	400cdm ²	
Contrast Ratio	1000:1	
Colour Bit Depth (LUT)	68.7 billion colors (12bit)	
View Angel	≥178°	
LUT	GAMMA2.0,GAMMA2.2,DICOM	
Video Input	DVI-I,VGA,CVBS,S-VIDEO,YPbPr,RGBS,SDI	
Video Output	DVI-I,CVBS,S-VIDEO,YPbPr,RGBS,SDI	
Optional	SDI	
Power Output	DC5V1A (USB-A)	
Power Requirements	100~240V, 50/60Hz	
Max Power Consumption	96W	
Typical Power Consumption	≤ 70W	
Cabinet Color	WHITE	
Dimensions	565.3x460.3x245mm	
Dimensions (Without Stand)	565.3x375.3x65mm	
Net Weight	6.5kg	
Net Weight (Without Stand)	9kg	
Hole Spacing	VESA: 100*100mm	
Certifications	CCC,CE	
Osd languages	English.	







ENDOTrolley

Mobile trolley for Endoscopy Systems

EndoTrolley is a versatile endoscopy system trolley. The ergonomically designed trolley offers plenty of space for your entire endoscopy suite and makes it easily accessible for a smooth workflow during procedures. The sturdy wheels make it easy to move the trolley when required.

99

Features

- Ergonomic design.
- · Vertically adjustable trays.
- Special handles for safe & easy movement.
- Secure hanger for enodscopes.
- Sturdy frame for Monitor placement.
- Compact & Portable

Specifications:

	Dimensions	1430 mm (H) x 510 mm (D) x 440 mm (W)
	Weight	40 kg
	Maximum Load Capacity	Top tray: 15 kg
		Per shelf: 20 kg
		Base Panel: 30 kg

Castors:

100 mm diameter
2x braked conductive twin-wheeled castors
2x braked twin-wheeled castors

Construction:

Workstation	Hoop: Mild steel		
	Shelves: Mild steel		
	Handles: Stainless steel 304		
	Side hanging bar: Stainless steel 304		
	Endoscope hanger: Acrylic		













Top Features



User-friendly



Full HD Image Viewing



Image Capture



View images in large size



Movie recording, editing and merging



Video Capture



ENDOClinic

Advanced Medical Image & Record Management System



Easy patient record search



Secure database



Gastrointestinal Endoscopy



Pulmonary Endoscopy



Urology Endoscopy



Laryngoscopy



Import/Export Image & Video



Backup patient data in CD/DVD/Pendrive



Laparoscopy



Hysteroscopy in IVF Centers



Veterinary Endoscopy



Create & edit new templates



Patient data summary in excel





Index




ENDOClinic LAN

Collaborate seamlessly with LAN Connectivity



105

ENDOClinic Connectivity



106







ENDOClinic Reports

ENDOClinic trusted by doctors for well-formatted customised reports

Sample Reports









ENDOClinic Feature

Increase Productivity with **ENDOClinic**

Features	EndoClinic 8iHD	EndoClinic 7	EndoClinic 8iHD LAN	EndocCinic 7 LAN
LAN Connectivity	X	×	✓	✓
Full HD image & video quality	✓	×	✓	×
Multi-user environment	X	×	✓	✓
Export/Import images & videos	✓	✓	✓	✓
Well-formatted multi-pages report	✓	✓	✓	✓
Patient data summary in excel	✓	✓	✓	✓
Easy patient record search	✓	✓	✓	✓
Secure database	✓	✓	✓	✓
Intelligent image cropping	✓	×	✓	×
Image & Video capture resolution	1920/1080p	720/576p	1920/1080p	720/576p
Compatibility	Windows 10 onwards	Windows 10 onwards	Windows 10 onwards	Windows 10 onwards
Required configuration	i5/i7 Processor, 8 GB RAM, 1 TB Hard disk, USB/PCI slot	i3/i5/i7Processor, 8 GB RAM, 500 GB Hard disk, USB/PCI slot	i5/i7 Processor, 8 GB RAM, 1 TB Hard disk, USB/PCI slot	i3/i5/i7Processor, 8 GB RAM, 500 GB Hard disk USB/PCI slot

High Definition

Standard Definition


ENDOClinic 8iHD LAN

ENDOClinic 7
ENDOClinic 7 LAN





SELF-EXPANDING METALLIC STENTS



Esophageal Stents
Biliary Stents
Pyloric Stents

Tracheobronchial Stents
Pancreatic Stents
Colo-rectal Stents

Urethral Stents

BLOOD BAGS & BLOOD BANKING EQUIPMENTS



LABSIDE Filter Bag



Inline Filter Bag



Triple Blood Bag







Blood Collection Monitor mini (m.bcm)



Double Pan Balance (m.DP)

Index



109

110





PERITONEAL DIALYSIS



APD m.Cycler



PD Twin Bag



Tri-Choice



4GYP



111

Our Product Lines

- Leukofilter Bags
- Blood Bags
- Blood Banking Equipment

Peritoneal Dialysis Systems

Non-Cardiac Metallic Stents

Flexible IV Bags Solutions

Video Endoscopy Systems

Endoscopy Ancillaries & Accessories

Endoclinic

www.mitraindustries.com

www.ottomed.com

Certificates





Mitra Industries Pvt. Ltd.

14/4, Delhi-Mathura Road, Faridabad, Haryana-121 003. INDIA

Telephone: +91-129-2271069 / 4099300

Toll Free: 1800 1020730

Email us @: customersupport@mipl.in Visit our Website: www.mitraindustries.com



CERTIFIED



Mitra Industries Pvt. Ltd.

B-226, Okhla Industrial Area, Phase-1, New Delhi-110 020. INDIA Toll Free: 1800 1020 307



Email us @: ottomed@ottomed.com, Visit our Website: www.ottomed.com

Ottomed **ENDOSCOPY**



www.ottomed.com









Manufactured by:

MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh -173 220. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com