


# mBLU ATLAS





**mBLU1**



**mBLU2**



**White Light**



**White Light**



**Dr. Meenakshi, Dr. Vinit Kumar and Dr. Prasad Bane**  
LTMG Hospital, Mumbai GI Surgery Department

**Procedure:** Video OGD Scopy

**Patient Gender:** Male

**Patient Age:** 49 years

**Symptoms:** Epigastric pain and regurgitation.



### **Pedunculated polyp at pylorus**

49 years male with epigastric pain and regurgitation with no significant findings on clinical examination unresponsive to empirical therapy for 15 days undergone diagnostic OGD scopy which showed pedunculated polyp at pylorus, biopsy taken was suggestive of hyperplastic polyp.

Follow up endoscopy was done after histopathology report to excise the polyp with snare and cautery. 1:10000 adrenaline in saline injected at the base and snare was used with cautery at the pedicle for excision. Polyp was retrieved using basket.

Post procedure there was no bleeding from mucosa at the base of excised polyp.


**mBLU**  
Case Studies



**White  
Light**



**mBLU 1**



**mBLU 2**

**Dr. Nimesh Shah, M.S., Mumbai**

(Consultant Laparoscopic & GI Surgeon & Endoscopist)

**Procedure:** Video Colonoscopy

**Patient Gender:** Male

**Patient Age:** 65 years

**Symptoms:** Bleeding in stool

**Impression:** Normal Colonoscopy

**Rare finding of long based Lipoma**

65 years old male presented with acute GI bleed with frank blood in stools. Patient was on Ecospirin 150mg till 2 days prior to admission. A colonoscopy & upper GI scopy was performed. Colonoscopy showed clots & streaks of blood in the colon. Upper GI scopy showed a long-based lipoma (confirmed by CT scan) which gave a double lumen appearance in the first part of duodenum which is a very rare observation. A CT angiography did not elicit active blush in any part of the bowel.

The patient settled with blood transfusion bowel wash & conservative management.

**mBLU**  
Case Studies



mBLU 1

## **mBLU** Case Studies

**Dr. A. S. A. Jegannathan, M.D., Madurai Scopy**

**Procedure** : Video Gastroscopy

**Impression** : Fundal Varices



White Light



mBLU 1

## mBLU Case Studies

**Dr. A. S. A. Jegannathan, M.D., Madurai Scopy**

**Procedure** : Video Sigmoidoscopy

**Patient Gender** : Male

**Patient Age** : 38 years

**Impression** : Plaque like lesions seen in sigmoid colon

## 2 Tech-in-1



### What does *mBLU* technology do?

Deep detection of mucosal and sub-mucosal tissue helps in better diagnosis and impacts the patient treatment outcomes. Option of comparing different tissue structures impacts diagnosis.

Introducing the




Experience

**mBLU's** advanced optical & light sensor technology used multi-LEDs (wavelength 400-700nm & 450nm) eliminating energy consumption & expensive Xenon light sources.

The possibility of detecting and characterising lesions / polyps increases due to brighter & clearer endoscopic images provided by **Image-Enhanced-Endoscopy (IEE)** using **mBLU**.

**mBLU** is MITRA's next-generation image technology which helps medical experts detect early-stage cancer and provide high quality **lifecare**.




## **mBLU1** Spectral Mucosal Tissue Structure Enhancement

Spectral Imaging assists in increasing the diagnostic yield and provides a better characterization of lesions.

mBLU1 provides an option of detection of lesions of “minimal” Esophagitis, of Dysplasia in Barrett’s Esophagus and of Squamous Cell Esophageal Cancer.

Being less invasive and time consuming than the classic dye methods, the use of mBLU1 aims to result in an improvement in the visualization of the early gastric cancer and improve characterization of colonic lesions.



## **mBLU2**

### Vascular Architecture Visualisation

mBLU2 Blue Light technology is vital to visualise the vascular microarchitecture on the mucosal surface. The blue light penetrates the superficial epithelium & visualises microstructures in the mucosa as well as the sub-mucosa. This mode can facilitate differentiation of adenoma from invasive cancer. The mBLU2 is dedicated to more close-up view or detailed view.



## **mBLU2 Mode**

mBLU2 is based on the principle of optical spectral visualisation using LED-At-Tip as light source eliminating the need of a Xenon light source.

mBLU2 does not need mechanical and color filters for generating spectral light as they have a tendency to deteriorate over a period of time. Deteriorated fiber bundles in scopes further reduce light intensity causing further darkness in the scope.

mBLU2 works with a single custom-made Fibreless LED-At-Tip with wavelength which is an average of blue and green filters. mBLU2 highlights capillary vessels of the subepithelial layer of the mucosa and provides better visibility of vessels of submucosal layer. We call it the “Depth” mode. (Use near focus mode while using mBLU2). mBLU2 is based on principle of absorption of particular wavelength of light by the hemoglobin.

Experienced endoscopists can derive good diagnostic results using mBLU2.



### **Advantages:**

- Precision diagnosis & identification helps in early detection of Cancer & Mucosal pattern changes.
- Both **mBLU** modes options are available in a single endoscopy unit with one touch button operation.
- Low operational cost & user-friendly.
- The diagnostic ability using the **mBLU** system allows more accurate discrimination of adenomas from non-adenomatous lesions as compared with conventional white light.

**mBLU1: Spectral Mucosal Tissue Structure Visualisation**

**mBLU2 : Deep Vascular Architecture Visualisation**



ON HOME	mBLU1 mBLU	ON LIGHT	IN PROCESS WHT BAL
MEDIUM ENHANCE	MODE 1 EXT IN	LOW PUMP	PEAK IRIS

Smart Touch Screen Interface

# OTTOMED Endoscopy Systems with mBLU

**Good Design provides a Great User Experience.**

**Simplified User Interface:** User Friendly operations via a single Smart Simple TouchScreen



MITRA is the first & the **only** Indian company to completely design and manufacture Endoscopes and related ancillary products in India in a state-of-the-art manufacturing facility.

Launched in 2012, Ottomed Endoscopes have become leaders in endoscopy technology having introduced the world's first endoscopes "LED-AT-TIP" patented technology.

Today, Mitra has successfully installed more than 1,000 endoscopy devices across the country.

Our specialised customer support teams are ready support your needs of prompt & friendly service.

**Thank you**

Mitra Industries Pvt. Ltd.

B-226, Okhla Industrial Area,  
Phase-I, New Delhi 110020.

Toll Free: 1800 1020 307

Email: [ottomed@ottomed.com](mailto:ottomed@ottomed.com)