

Insight^{Plus}

HD Video Processor



Model : OVP-IN Plus



COMPATIBLE WITH

- Bronchoscope
- Pleurascope
- Laryngoscope
- Ureterscope
- Cystoscope



Video Recording



Image Capturing



Compact &
Light Weight

Features:

- **In-Built Image Capturing & Video Recording System**

In-built Image & Video recording system stores captured Images & Videos directly to an External pen drive

- **5-in-1 Connectivity - Broncho/Pleura/Laryngo/Uro/Cysto Scopes**

5-in-1 video processor compatible with Broncho/Pleura/Laryngo/Uro/Cysto endoscope models.

- **Freeze Function**

Image freeze function helps select clearest Image.

- **Single Click Connector**

Single-handed Single-Click Connector makes the Endoscope more reliable & user-Friendly.

- **Compact & Light Weight**

LED-At-Tip technology helps make the system Compact & Light weight. Easy to carry anywhere.

Specifications:

Model	OVP-IN Plus	
Power Supply	Voltage	100-240 VAC
	Voltage fluctuation	Within±10%
	Frequency	50/60 Hz
	Frequency fluctuation	Within±3 Hz
	Consumption electric power	17VA
	Fuse rating	3 A, 250 V
	Fuse size	ø5.0 x 20.0 mm
Size	Dimensions	300(W) x 295(L) x 80(H) mm
	Weight	4.8 kg
Classification (medical electrical equipment)	Type of protection against electric shock	Class 1
	Degree of protection against electric shock of applied part	Type BF applied part. Where no classification mark appears, the device is a TYPE BF applied part.
	Degree of protection against explosion	The Video Processor should be kept away from amammable gases.
Observation	Video signal output	DVI-D x 2, HD-SDI
	White balance	By pressing the white balance switch on the front panel, automatic white balance can be performed.
	Zoom	1.2X & 1.5X
	Image enhancement mode setting	Fine patterns or edges in the endoscopic images can be enhance electrically to increase the image sharpness.
Portable Memory	Media	USB mass storage
	Image recording	JPEG format
	Movie recording	MP4
Illumination	LED (White)	At Endoscope Tip

Manufactured by:

MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan,

Himachal Pradesh - 173 220. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

