



Qubey-II

Fully Automatic Twin Scope Endowasher System

360° Auto Disc Rotation Water Jet




Minimize your Hygiene & Safety concerns. Disinfect your scope in the **Qubey-II Endowasher** after every procedure to maximize safety of the Patient, Staff and yourself.











Features

- Fully automatic operation.
- Touch Screen user friendly Interface for ease of operation. 2.
- Memory for 50 cycle storage. 3.
- Inbuilt printer for washing and disinfection report. 4.
- USB port for data retrieval.
- 6. Leakage testing facility.
- Visual indicator for disinfection chemical level. 7.
- Specially designed tub to ensure minimum wastage of disinfectant & detergents. 8.
- Powerful channel drying mechanism post cleaning & disinfectant cycle.
- 10. Can be operated in manual mode as well.
- 11. Compatible with all Endoscopy systems.

Specifications - Ottomed Fully Automatic Endowasher (OEW-2)

Ambient temperature		y Automatic Endowasher (OEW-2) 10°c - 40°c
Relative Humidity		30% to 85%
Water Supply Requirement	Flow Rate	12L/min – 16L/min
	Pressure	0.2MPa - 0.5 MPa
Main Unit Specification		
No. of Accommodated Endoscopes		2 nos.
Disinfectant Tank Capacity		17L
Detergent Tank Capacity		1.5L
Leakage Testing		Bubble detecting during immersion (Visual)
Operating Modes		
Auto Mode		
ν 	Operations	Time Range
	Wash cycle	1 – 10 min
	Disinfectant time	1 – 99 min
	Dry time	1 – 10 min
	Detergent wash time	1 – 10 min
Alarm		At Low level
		End of process
Display		Touch Screen Panel
Dimension		600 x 600 x 600mm (L x W x H)
Weight		60 Kg Approx
Power Supply		
Input Voltage		220 V
Allowable Voltage Fluctuation		±10%
Input Current		1.5 A
Frequency		50Hz
Allowable Frequency Fluctuation		±3









MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com