



# Video Pleuroscope

Insight plus Series

A Slim Size Thoracoscope with Fiberless LED at Tip technology




Model: OTS-IN+



LED At Tip



CMOS sensor




Laser compatible



Single click connector


**Features:**

- LED-At-Tip Technology for longer scope life.
- CMOS Sensor provides unmatched image quality to explore smallest details.
- Ergonomic, very slim insertion tube diameter 6mm helps smooth insertion & improves procedures performance.
- Flexible angulation, UP 180° & DOWN 130° enables easy orientation inside the pleural cavity.
- Light weight design for easy handling.
- Single Click Connector.



**Specifications:**


Model	OTS-IN+	
Optical System	Field of view	120°
	Direction of view	Forward View
	Depth of field	3-60mm
Insertion Tube	Distal End outer diameter	6 mm
	Rigid Section outer diameter	6 mm
	Maximum Insertion section outer diameter	6.7 mm
	Distal End enlarged	1. Illumination LEDs (W) 2. Image Sensor 3. Instrument Channel Outlet
Instrument Channel	Working Length	300 mm
	Channel inner diameter	2.8 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endotherapy accessories enter and exit the endoscopic image	
Bending Section	Angulation range	UP 180°, DOWN 130°
Total Length		595 mm
Laser Compatibility		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-IN Plus



Distal End Tip



Light Weight Design



Flexible Angulation



Manufactured by:  
**MITRA MEDICAL SERVICES LLP**  
 Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA  
 Toll Free No.: 1800 1020 307  
 E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:  
**MITRA MEDICAL SERVICES LLP**  
 B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA  
 Toll Free No.: 1800 1020 307  
 E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

\*Specifications are subject to change.

Printed - 03/JAN/2022