



ENDOPLUS HD

Full HD Laproscopic Camera System



Model: OVP-EP



Clear, Crisp & Brighter Imaging













Zoor



Remote Switch



Compact & Light Weight



Single Click Connector

• Full HD Imaging

High Definition CMOS Sensor gives Clear, Crisp & Bright imaging with outstanding flicker free quality.

• Inbuilt image capturing

In-Built Image capturing system captures images directly to an External pen drive.

Zoom Function

Image zooming with 2 steps 1.2x & 1.5x, improves visualisation during Minimally Invasive Surgery (MIS).

• Compact & Light Weight

Compact size & Light weight system makes it easy to move & carry anywhere.

Remote Switch

Assigned Remote switches on camera head for easy operation.

• Easy scope attachments

Adjustable camera head can connect to most of the branded rigid & flexible endoscopes.

• Single Click Connector

Single handed Single Click Connector makes the camera more reliable & user-friendly.

Specification:

Model	OVP-EP	
Operating Environment	Ambient Temperature	10°c - 40°c
	Relative Humidity	30% to 85%
Power Supply	Voltage	100V to 240V
	Voltage Fluctuation	±10%
	Frequency	50/60Hz
	Frequency Fluctuation	±3Hz
	Current Rating	0.7A-0.3A
	Fuse Rating	4A, 250V
TV system	Image pick-up element	CMOS- 2MP
Picture quality	Resolution	Horizontal 1080 lines
	S/N	49 dB min.
	Compatibility	Camera head is compatible with Ottomed HD Video System. Auto white balance should be adjusted.
Safety	Classification of medical electronic	Protection level of the fitting part BF type.
	apparatus	
Mounting	Mounting method	Screw mounting M22x0.75
Optics	Elimination of moiré	HD video camera head is available with built-in moiré elimination filter.
Waterproof	Level	Watertight when fitting the video adapter to the HD video camera head and fitting the waterproof cap to the camera connector.
Resistance to disinfection	Cleaning, Disinfection	Water disinfect by fitting the video adapter to the HD video camera head and fitting the waterproof cap to the camera connector.
	Remote Switch 1	Zoom
Remote Switch	Remote Switch 2	Record
	Remote Switch 3	White Balance
Observation	Video signal output	HDMI, HD-SDI &VBS, Simultaneous output is possible.
	Enhancement	Fine patterns or edges in the endoscopic images can be enhanced electrically to increase the image sharpness.
	Switch enhancement mode	The enhancement level can be selected from 3 levels (Low, Med and High) using the image enhancement mode on the front button panel.
Storage	External Storage	USB Memory Stick
Size	Dimensions	305(W) x 235(L) x 80(H) mm
	Size	3.5kg approx.









Manufactured by:

MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Printed - 21 JULY 2022