



# ULTRA-SLIM VIDEO GASTROSCOPE S2+ Series







LED At Tip

## Model: OVG-S2+US



CMOS Sensor

**Remote Button** 



Single Step Water Proof Connector

### Features:

- Super Bright LED-At-Tip: Experience 30% more brightness & 15% wider spread, thanks to the new Phosphor coating on the LED substrate, providing high reliability over years of use.
- CMOS Sensor provides unmatched image quality to explore smallest details.
- Ultra-Slim insertion tube unit.
- Remote buttons for Easy operation
- Ultra-light weight, Single Step Waterproof Connector makes the Endoscope more reliable & user-friendly.









### **Specifications:**

Model	OVG-S2+US	
	Field of view	120°
Optical System	Direction of view	Forward View
	Depth of field	3-100mm
	Distal end outer diameter	7.8 mm
	Insertion tube outer diameter	7.8 mm
Insertion Tube	Distal end enlarged	Air/water nozzle Up Left Objects lens Down
	Working Length	1055 mm
Instrument Channel	Channel inner diameter	2.2 mm
	Minimum Visible Distance	3.5 mm from the distal end
	Direction from which endotherapy accessories enter and exit the endoscopic image	
Bending Section	Angulation range	UP 210°, DOWN 90°
		RIGHT 100°, LEFT 100°
Total Length		1365 mm
Compatible Video Processor		OVP-S2+



#### Manufactured by:

MITRA MEDICAL SERVICES LLP Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

#### Distributed by: MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA Toll Free No.: 1800 1020 307 E-Mail: ottomed@ottomed.com, Website: www.ottomed.com Printed - 21 JULY 2022