




ULTRA-SLIM VIDEO VETERINARY ENDOSCOPE S2+ Series



Model: OVV-S2+US



LED At Tip



CMOS sensor




Remote Button



Single Step,
Water Proof
Connector

Features:

- Multi-beneficial super bright LED-AT-TIP Technology for longer scope life.
- CMOS Sensor provides unmatched image quality to explore smallest details.
- Easy manoeuverability.
- Remote buttons for easy operation.
- Ultra light weight, single step waterproof connector makes endoscope more reliable & User-friendly.



Specifications:

Model	OVV-S2+US (I/S)	
Optical System	Field of view	120°
	Direction of view	Forward View
	Depth of field	3-100mm
Insertion Tube	Distal end outer diameter	7.8 mm
	Insertion tube outer diameter	7.8 mm
	Distal end enlarged	
Instrument Channel	Working Length	I-1400 mm/ S-1055 mm
	Channel inner diameter	2.8 mm
	Minimum Visible Distance	3 mm from the distal end
Bending Section	Direction from which endotherapy accessories enter and exit the endoscopic image	<input type="checkbox"/>
		<input checked="" type="checkbox"/>
Total Length	Angulation range	UP 180°, DOWN 120°
		RIGHT 120°, LEFT 120°
Compatible Video Processor		OVP-S2+

*Specifications, Designs and Models are subject to change without any obligation on the part of the manufacturer



Manufactured by:

MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA

Toll Free No.: 1800 1020 307

E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Printed - 18 FEBRUARY 2022