



VIDEO URETEROSCOPE

Insight Plus Series





Model:OUR-IN+



LED At Tip



CMOS Sensor



Laser Compatible



Single Click Connector

Printed - 06 NOVEMBER 2023

Features:

- Super Bright LED-At-Tip: Experience 30% more brightness & 15% wider spread, thanks to the new Phosphor coating on the LED substrate, providing high reliability over years of use.
- CMOS Sensor provides unmatched image quality to explore smallest details.
- Ergonomic, Slim and Flexible for Smooth Insertion.
- Easy Manoeuverability.
- Laser Compatible.
- Ideal for Lithotripsy.
- Light weight for easy handling.
- Flexible angulation (270° Up & Down) enables optimal visualization in the Upper & Lower Calyxes respectively.
- Single Click Connector.









Specifications:

Model		OUR-IN+
Optical System	Field of view	120°
	Direction of view	Forward View
	Depth of field	3-60mm
Insertion Tube	Distal end outer diameter	2.85 mm Bullet Shape
	Insertion tube outer diameter	2.85 mm
	Distal end enlarged	1. Illumination LEDs (W)
		2. Image Sensor
		3. Instrument Channel Outlet
	Working Length	670mm
Instrument Channel	Channel inner diameter	1.2 mm
	Minimum Visible Distance	3 mm from the distal end
	Direction from which endotherapy accessories enter and exit the endoscopic image	
Bending Section	Angulation range	UP 270°, DOWN 270°
Total Length		965 mm
Laser Compatibility		Nd: YAG, 810nm Diode
Compatible Video Processor		OVP-IN2



Manufactured in India by:

MITRA MEDICAL SERVICES LLP Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA Toll Free No.: 1800 1020 307 E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by: MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA Toll Free No.: 1800 1020 307 E-Mail: ottomed@ottomed.com, Website: www.ottomed.com