

ENDOPUMP-II

Greater Visibility for improved Diagnosis & Detection



RECOMMENDED APPLICATIONS:

- Auxiliary water jet/Biopsy channel irrigation for gastric washing, bleed etc.
- Auxiliary water jet irrigation for washing colonic mucosa.
- Quick filling for Endoscopic Ultrasound procedure.

**Peristaltic
Irrigation Pump**

Features:

- Reverse Flow Technology for precise Flow control
- In-built Safety Alarm.
- LED Indications.
- Noiseless operation.
- Easy selection of water flow rate.
- Foot switch control.
- Universal Compatibility.



Specifications: m.EP-II

Operating Environment

Ambient Temperature	10°C – 40°C
Relative Humidity	30% to 85%

Main Unit Specifications

Power Supply

Input Voltage	100V - 240V
Allowable Voltage Fluctuation	±10%
Current	0.4A – 0.2A
Frequency	50Hz / 60Hz
Water Flow Range	100ml/min – 900ml/min
Foot Switch Type	Pedal Type
Flow Indicator	LED

Dimensions

Length	255mm
Width	265mm
Height	150mm
Weight	Approx. 2.7 Kg Without Water Tank



Made in India by:

MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan,
Himachal Pradesh - 173 220. INDIA
Toll Free No.: 1800 1020 307
E-Mail: ottomed@ottomed.com, Website: www.ottomed.com

Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA
Toll Free No.: 1800 1020 307
E-Mail: ottomed@ottomed.com, Website: www.ottomed.com