




ENDOCARBON-III Co₂ Insufflator



Efficient & Compact machine for Patient comfort.

Advance micro-controller technology. Special customized & pre-programmed single touch mode for POEM procedure.

*Actual Gas Input Pressure Indicator *Light Weight *Compact Size *Universal compatibility *Single touch flow selection




Main Features

- Cr Time Saving and Patient Friendly Examination
- Carbon Dioxide is absorbed by human tissue about 150 times faster than the conventional air.
- CO₂ Rapid Absorption properties keep abdominal distension and related pain to a minimum during the process and speeds up recovery afterwards.



Specifications- m.EC-III		
Operating Environment		
Ambient Temperature		10°c – 40°c
Relative Humidity		30% to 85%
Main Unit Speci cat	ions	
Pneumatics		
Gas		CO_2 (Medical Grade)
Input Pressure		1 bar - 7 bar
Flow Rate	POEM	1.0 L/min ±10%
	LOW	3.5 L/min ±10%
	HIGH	5.0 L/min ±10%
Power Supply		
Input Voltage		100V - 240V
Allowable Voltage Fluctuation		±10%
Power Rating		24VA
Frequency		50Hz / 60Hz
Allowable Frequency Fluctuation		±3 Hz
Dimensions		
Length		155mm
Width		255mm
Height		150mm
Weight		Approx. 2.3 Kg

System Chart





Made in India by: MITRA MEDICAL SERVICES LLP

Plot No.: 21-B, Sector-1, Parwanoo, Dist.- Solan, Himachal Pradesh - 173 220. INDIA Toll Free No.: 1800 1020 307 E-Mail: ottomed@ottomed.com, Website: www.ottomed.com Distributed by:

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1, New Delhi - 110 020. INDIA Toll Free No.: 1800 1020 307 E-Mail: ottomed@ottomed.com, Website: www.ottomed.com